# AILORIA

## LUSTRE

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Protective cap

### **Product information**

On/off switch with light indicator (charge indicator)

360° rotation

2 speed levels

Integrated suction function: all skin particles are automatically sucked in

Pleasant soft surface in diamond cut design

Charging time: 1,5 hours

Working time (when fully charged) : approx. 30 minutes

### **Charging**

Connect the USB charging cable to a USB interface (e.g. notebook, PC) or to the supplied USB charging adapter and connect it to a power source.

As soon as the cornea remover is fully charged, the charging indicator switches from red to green. Then disconnect the cornea remover LUSTRE from the power source, as overcharging may reduce the life of the battery.

## **Application**

Place the desired sanding attachment - coarse grit for thick cornea, fine grit for a soft finish - on the top of the product. Place it in place and press until it clicks into place.

Switch on the charged electric cornea remover by pressing and holding the on/off switch on the bottom for 2 seconds. The product has 2 speed settings: Press and hold the on/off button for 2 seconds for smooth rotation. Pressing the on/off button again will activate the higher speed setting.

Apply LUSTRE to the desired and previously cleansed areas of skin and remove calluses. Do not hold the cornea remover in the same place for more than 10 seconds. Do not apply strong pressure - thanks to the advanced technology, it only needs to be gently applied to the skin.

Once you have achieved the desired result, switch off the cornea remover with the on/off button on the bottom of the device by pressing and holding it again for 2 seconds.

## **Changing the sanding attachment**

Carefully pull the sanding head upwards with two fingers on the notches and remove it.

Place the desired sanding head in the centre and press it down until it locks into place.

## **Cleaning**

Carefully remove the grinding head from the unit. Clean it under warm water. Empty the container. The cornea remover must not be cleaned under running water. Thanks to its soft surface, LUSTRE can be easily cleaned with a soft, damp cloth.

## **Storage**

After use, store the cornea remover in a dry and clean place and protect it from sunlight.

## **Technical data**

Model: LUSTRE

Material: ABS

Dimensions: 6,6 x 6,6 x 9,8 cm

Weight: 134,8 g

Battery: 400 mAh Li-ION

Voltage: 3.7 V

## **Safety Instructions**

Please read the following warnings before using the product. Improper use can lead to injuries.

This appliance may be used by children over 8 years of age and by persons with reduced physical, sensory or mental abilities or lack of experience and/or knowledge if they are supervised or instructed in the safe use of the appliance and have understood the resulting dangers.

Please keep the product away from children.

The cornea remover is not waterproof and should not be used in any situation (e.g. shower) where it may come into contact with water. If water gets into the product, please switch it off immediately and stop using it.

The cornea remover is only suitable for corneal removal. To avoid injury, do not use the device on smooth, non-cornified skin.

Do not use the product on irritated, inflamed or infected skin or

open wounds as this may cause irritation or cuts.

If you are a diabetic or suffer from circulation problems, please consult a doctor before using the product.

The cornea remover must not be used immediately after an operation or on scars.

If you feel pain while using the product, discontinue use immediately.

Do not introduce foreign bodies into the product to avoid damage.

This product must not be disassembled under any circumstances.

## **Disposal**

According to the EU directive, this product must not be disposed of together with other household waste. Find out about and follow local regulations on the separate collection of electrical and electronic equipment. Proper disposal of old products helps to prevent negative effects on the environment and health.

## **Manufacturer**

VEHNS GROUP GmbH  
Theatinerstrasse 40-42  
80333 Munich  
Germany

[www.vehnsgroup.com](http://www.vehnsgroup.com)  
[service@vehnsgroup.com](mailto:service@vehnsgroup.com)  
[www.ailoria.de](http://www.ailoria.de)

MADE IN CHINA

Changes and errors excepted  
The manufacturer accepts no liability for damage caused by incorrect, improper or incompatible use of the product.

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© VEHNS GROUP GmbH
VEHNS GROUP GmbH, Theatinerstrasse 40-42