
DEPLAQUED

DE Bedienungsanleitung **EN** User manual **FR** Manuel de l'utilisateur **IT** Manuale di utilizzo **ES** Manual de Usuario **NL** Gebruiksaanwijzing **PL** Instrukcja obsługi **RO** Manual de utilizare **PT** Manual do utilizador **SV** Bruksanvisning **DA** Brugervejledning **CS** Uživatelská příručka **NO** Bruksanvisning **SK** Používateľská príručka **FI** Käyttöohje **LT** Naudotojo vadovas **BG** Ръководство за потребителя **HR** Korisnički priručnik

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
PRODUCT FEATURES

The AILORIA DEPLAQUED sonic dental cleaner is ideal for cleaning teeth at home. It removes tartar, plaque and discolouration, even from hard-to-reach areas.

- Intelligent high-frequency vibration
- Ergonomic design for comfortable handling
- 5 intensity levels
- Cordless battery operation
- Charging time approx. 2 hours
- Battery life approx. 6-8 hours when fully charged
- Charging of the lithium-ion battery via USB

COMMISSIONING

Charge the sonic dental cleaner before the first use. Charging time is approximately 2 hours to charge the sonic dental cleaner before first use or fully discharged battery.

CHARGING PROCEDURE

The connector for the charging cable is located at the bottom of the back of the unit. Open the silicone cap and connect the charging cable to the sonic dental cleaner. Connect the other end to a power source (wall socket) using the USB charging adapter provided. Alternatively, you can also connect the USB cable to a USB interface (e.g. PC or notebook).

The sonic dental cleaner cannot be used during the charging process. A red light shines in the transparent base for the duration of the charging process. When the charging process is completed, this is indicated by a green light.

APPLICATION

Press the ON/OFF button to switch on the sonic dental cleaner. Select the desired intensity using the arrow keys.

The following 5 intensity levels are available:

Intensity level 1: 8321 vibrations per minute

Intensity level 2: 9787 vibrations per minute

Intensity level 3: 11195 vibrations per minute

Intensity level 4: 12423 vibrations per minute

Intensity level 5: 13320 vibrations per minute

Start cleaning the teeth by gently placing the metal attachment on the areas with plaque and discolouration. Now slowly move the metal tip over the areas to be cleaned. This process can be repeated several times until the desired result is achieved. The sonic dental cleaner can also be used on the edges of the teeth and in the spaces between the teeth. Move the sonic dental cleaner carefully to prevent possible injuries to the gums. Repeat the application 1 to 2 times per week. The DEPLAQUED can also be used on dental crowns and bridges. Proceed carefully.

Clean the metal head with water after each use and wipe the housing with a damp cloth. Sterilize the sonic dental cleaner with alcohol if necessary.

SAFETY INSTRUCTIONS

- Always keep the USB cable and USB adapter away from water and other liquids.
- Use only the USB cable supplied with this product. Using other cables may damage the device or affect charging performance. Do not bend or kink the USB cable excessively. Make sure the power

adapter used with the product's USB cable is suitable and compatible.

- Stop using the product immediately if it malfunctions or if you notice any signs of a possible fire hazard.
- Disconnect the USB charging cable before cleaning the sonic dental cleaner.
- Do not use the product if it is damaged. The sonic dental cleaner may only be repaired by the manufacturer.
- Children under the age of 14 and persons with reduced physical or mental abilities, or with a lack of experience and knowledge, may use the product only under adult supervision and instruction. Always keep DEPLAQUED out of the reach of children.
- DEPLAQUED is intended solely for cleaning teeth and must not be used for any other purpose.
- If excessive gum bleeding or pain occurs during use and persists, consult your dentist.
- The product should not interfere with pacemakers. As a maximum precaution, however, we recommend consulting your doctor before use.
- If you have had jaw or gum surgery within the last two months, or if you have any medical concerns, consult your doctor before using DEPLAQUED.
- Clean the metal attachment with water after each use and wipe the housing with a damp cloth. If necessary, sterilize the dental cleaner with alcohol.
- The sonic dental cleaner is protected against splashing water according to the IPX6 standard.

However, do not hold the product under water for an extended period.

- Do not overcharge the battery. Disconnect the device from the charger once it is fully charged.
- Store the device in a cool, dry place and avoid extreme temperatures.
- Do not expose the battery to direct sunlight and do not store the device in damp environments.
- Protect the device from impacts, drops and strong pressure to prevent damage to the battery.
- Do not attempt to remove or replace the battery yourself.
- If you notice any signs of battery damage or leakage, stop using the device and contact customer service.
- Avoid contact with leaking fluid. In case of skin contact, rinse thoroughly with water immediately.
- At the end of its service life, dispose of the device in accordance with local regulations for electronic waste.
- Do not throw the battery into fire or water.
- If the device overheats, emits smoke or shows any other signs of a defect, switch it off immediately and disconnect it from the power source.
- Keep the device away from flammable materials until the issue has been resolved.

TECHNICAL DATA

Rated voltage: 3.7V

Rated current: 200mA

Battery capacity: 300mAh

Charging voltage: 5V


DISPOSAL INFORMATION

Waste electrical and electronic equipment must not be disposed of with household waste. You must therefore not simply throw away your old appliances. They must be recycled. In order to improve waste separation, the Electrical and Electronic Equipment Act therefore stipulates that all electrical and electronic equipment must be permanently marked with the corresponding symbol, the crossed-out dustbin. Free collection points for old electrical appliances are available for you to return them. You can obtain the addresses from your city or municipal administration. Batteries and accumulators that are not enclosed in the old appliance, as well as lamps that can be removed from the old appliance without being destroyed, must be removed from the old appliances before disposal. As a consumer, you are generally obliged to hand in used batteries and accumulators as well as lamps separately at the special collection points. Make sure that the batteries and accumulators are completely empty. Then carefully remove the batteries and accumulators from the old appliance.

The following batteries are included with this item:
1 rechargeable battery, type Li-Ion, 3.7V, 300mAh.

Subject to alterations and errors

The manufacturer accepts no liability for damage caused by incorrect, improper or incompatible use of the product.

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© VEHNS GROUP GmbH
DONG GUAN YI YUAN DA ENERGY TECHNOLOGY CO. , LTD
VEHNS GROUP GmbH, Theatinerstr. 40-42, 80333 Munich,
VEHNS GROUP GmbH, Theatinerstraße 40-42